

GO CHAIR



Small, Light and Portable.
The New revamped Go-Chair



*Shown in Ruby Red

The all-new go chair is re-engineered from the ground up, offering a sleek, bold new style available in an array of contemporary colours.

Enhanced performance and comfort, along with feather-touch disassembly, allows you to enjoy light-weight travel and independence on the go.

With an increased weight capacity of 21 stone 6 lbs and a sleek, bold look, the new Go-Chair makes travel easy.

Under seat storage



Feather-touch disassembly



Colour Options



Amber Orange



Rose Quartz



Sapphire Blue



Citrine Yellow



Pearl White

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —



Tel: 01869 324600
Email: sales@pride-mobility.co.uk
Web: www.pride-mobility.co.uk

Features



- Choice of 6 colours
- 45 cm x 42 cm, 5-point swivel seat with latch release
- Increased 21 stone 6lbs, weight capacity
- Faster speed - up to 3.7 mph
- Larger, 18 amp batteries
- Improved agility - longer wheel base; wider stance; compact turning radius of 63 cm (25.4")
- Dual swing away, under seat storage
- Larger, height & width adjustable armrests
- Larger footplate
- One hand feather-touch disassembly

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Seat to ground height	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
Stone								Stone	Stone		
21 Stone 6lbs (136 kg)	3.7 Mph (5.9 km/h)	4 cm (1.6 inch's)	64 cm (25.4 inch's)	84 cm (33 inch's)	57 cm (22.6 inch's)	Front: 15 cm Rear: 23 cm	56 cm - 58 cm	30.39 kg	16.32 kg (rear)	Folded Seat 13.6 k.g	x2 18 Ah Deep Cycle (14 kg each)

Warranty

Lifetime limited warranty on frame.
1-year limited warranty on electronics.
18-month warranty on drive control system.
2-year limited warranty on drivetrain.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.

#myprideproduct



@PrideMobilityUK

#socialmobility